



WoundStat™
HEMOSTATIC WOUND TREATMENT

Technology designed to save the lives of those who serve our country

- » Uncontrolled hemorrhage is the primary cause of death in battle
- » 80% of wounds suffered by warfighters today are severe, irregular wounds caused by IEDs
- » Until now, NO product has been proven to consistently stop high pressure bleeding...

Now, that product is here.
WoundStat stops bleeding
in seconds, not minutes.

“...20% of deaths are still regarded as preventable. The primary issue in these deaths is uncontrolled hemorrhage.”

J. Holcomb, 2006, Presentation at Advanced Technology Applications for Combat Casualty Care (ATACCC)

WoundStat™

HEMOSTATIC WOUND TREATMENT

- ★ **Stops** high pressure **bleeding in seconds, not minutes**
- ★ Four independent studies show **WoundStat is the ONLY hemostatic agent to demonstrate 100% survival**
- ★ U.S. Army testing shows WoundStat had the **least post treatment blood loss**
- ★ **Longest duration** of hemostasis
- ★ **No exothermic reaction**
- ★ Dense granules **resist wind and rain**
- ★ **Safe to use** – no transmission of bloodborne pathogens
- ★ **Easy to apply** by self or buddy/medic (Per TCCC)
- ★ **Ready to use** – no mixing
- ★ **Repack and remold existing product** if re-bleeding occurs
- ★ Demonstrated **stability for up to 3 years**
- ★ **Affordable**



APPLICATION INSTRUCTIONS



Tear open WoundStat pouch.



Empty contents directly into wound.



Firmly pack into all areas of the wound and hold pressure for up to 3 minutes. If rebleeding occurs, repack existing WoundStat.



Wrap and secure a bandage to maintain pressure during transport. Move wounded quickly.

**Rapid.
Ready.
Reliable.**

**TO LEARN MORE /
TO PLACE AN ORDER,
PLEASE CONTACT:**



XSPEC TECHNOLOGY SDN BHD

No.8, Jalan Industri USJ 1/8
Taman Perindustrian USJ 1
47600 Subang Jaya, Selangor
MALAYSIA

Tel: +603 8023 1161

Fax: +603 8023 1167

E-mail: inquiry@xspec.com.my

TraumaCure

For more information visit...

